

CP 550



EXTENSIVE RANGE OF TESTS TO CARRY OUT COMPREHENSIVE EYE EXAM IN A SIMPLE WAY!



remote
control



CP 550

CHART PROJECTOR



Min : 2 m



Max : 7 m



pola

Optometry





EXTENSIVE CHOICE OF TESTS

- Visual acuity: Letter, Numbers, Snellen, Landolt C, children charts.
- Binocular balance and special tests: Astigmatism, polarised tests, red/green tests, Aniseikonia, Schober, Phoria...
- 2 types of progression available : A type or B type.

FLUID ACUITY ASSESSEMENT

- Can be managed either by remote control or by APH550 automatic phoropter keyboard.
- Long lasting LED illumination provides bright and clear images.
- Fast and quiet disk rotation ensures patient comfort.

OPTIONS

- Wallmount bracket, telescopic tube, polarised filter glass, R&G glass

PROJECTION	TESTS (A OR B TYPE)
Size	335 x 222 mm (at 5m)
Magnification	30 x (at 5m)
Distance	Projection distance 2.m~7m
NUMBER OF CHARTS	50
PROJECTION SPEED	0.03 seconds
AUTOMATIC SHUT-OFF	5mn
PROGRAM SETTINGS	2 specific programmable test sequences
STANDARD ACCESSORIES	Table base stand Remote control + Batteries Spare fuses 250v - 2AT Projection screen Calibration sheet Test card Dust cover
	Optotypes
	Numbers, Letters, Landolt C, Snellen, Children
	Mask and tests
	Lines (Horizontal/Vertical) Astigmatism, Cross-cylinder, Jackson Cross, Worth 4 dots Binocular balance Duochrome balance Coincidence Schober, Phoria Minute stereo Fixation point Eye picture
	DIMENSIONS AND WEIGHT
	248(W) x 234(H) x 336(D)mm - 4.5kg
	PROJECTION LAMP
	30 VA
	POWER SUPPLY
	AC 100 - 240V, 50/60Hz
	CE Conformity marking

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.

PRODUCT FUNCTIONS



Illumination



Projection distance



Linear polarisation



Wireless connection capability