


CS 550



# HIGH DEFINITION SCREEN WITH RED AND GREEN CHARTS



CS 550


CHART SYSTEM

Min : 1,5 m  
Max : 6 m



Optometry





## HIGH QUALITY VISUAL ASSESSMENT

- Over 100 kinds of general and special charts
  - Visual acuity: Snellen, Landolt c, symbols for children, numbers and letters including ETDRS.
  - Special tests: contrasts, binocular vision, aniseikonia, red/green, astigmatism...
- Optotypes randomisation and reversible tests for mirrored room.
- 2 types of progression available : A type or B type.

## OPTIMISED REFRACTION CHAIN

- Can be managed either by remote control or by APH550 automatic phoropter keyboard.

## OPTIONS

- Floor stand, R&G glass

<b>SCREEN SIZE</b>	19"	<b>TEST CHARTS</b>	Type A (Decimal indication) Type B (Feet and Metric indication)
<b>LUMINOSITY</b>	80 – 320Cd/M <sup>2</sup>	<b>Optotypes</b>	Letters, Landolt C, Numbers, Letters, Snellen, Children charts,
<b>CONTRASTS</b>	From 100% to 1.5%	<b>Mask and tests</b>	Astigmatism Binocular vision Coincidence (aniseikonia) Red/Green Heterophoria Stereoscopic vision
<b>PROJECTION DISTANCE</b>	1.5 to 6m		
<b>NUMBER OF CHARTS</b>	100		
<b>PROGRAM SETTINGS</b>	2 specific programmable test sequences		
<b>STANDARD ACCESSORIES</b>	Mini USB Cable for recharging remote control, Wall mounting bracket	<b>DIMENSIONS AND WEIGHT</b>	436(W) x 362(H) x 60(D)mm - 5.5kg
		<b>POWER SUPPLY</b>	AV 100-120V/220-240V, 50-60Hz
		<b>CE</b> Conformity marking	

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.

### PRODUCT FUNCTIONS



Screen Size



Projection distance



Wireless connection