

A typical in-store installation

Essibox

Information link between the operator and Essibox system, communicating via discreet luminous displays.

- Weight: 1,2 kg • Size: 15x15x15 cm • Power supply (voltage): 220 V-50/60 Hz
- Electrical consumption: 15 W • Communication through radio waves • CE marking conformity

Essibox Server

Manages Essibox applications, communications and machines. Detailed specifications contained in the sales terms.

Guarantee and related services

Complete service including installation, training, support, equipment guarantee, anti-virus protection and users' rights for the applications. Detailed specifications contained in the sales terms.

Operating conditions

The computers present in your shop and on which the Essibox applications will run (Up to 5 simultaneous users), must be equipped with:

- The Microsoft® Windows 2000, Windows XP or Windows Vista operating system;
- Microsoft® Internet Explorer version 6 or 7.

The shop should also have

- A permanent, "always-on" ADSL link. Essibox architecture needs access to the web in order to operate.
- A local network (cable or Wi-Fi).

Meets the ISO 16284 standard [OMA 3.07 compatible].

To allow for ongoing improvements, these specifications are non-contractual and are subject to modification without notice.

ORGANISING
VISION



ESSILOR INTERNATIONAL Instruments Department

64 bis, avenue Aubert – 94306 Vincennes Cedex France - Tél. : 33 (0)1 72 94 71 00 – Fax : 33 (0)1 72 94 70 27
www.organising-vision.com

Essibox.
6 new ways to make your working day easier!






1. Essibox guarantees continuous data transfers

Each day, you collect and process information from a range of mechanical and computerised equipment (optometry, measurements, store management, frames, etc.). Although a number of channels already exist, the transfer of information between each of your items of equipment is often patchy and fragmented. This lack of automation separates each task, requiring you to enter data several times unnecessarily, and can be very time consuming.

With Essibox, Essilor offers you the possibility to connect up all of your tools, creating a single continuous flow of information. A genuine control centre in which the various components of your store are unified, Essibox centralises all data flows, making you more efficient and giving you more time to concentrate on your clients.



**FEWER
RISKS OF
ERRORS**

**SAVES
YOU TIME**



2. Essibox blends in seamlessly with your organisation

Requiring no modifications to your existing organisation, Essibox is directly incorporated within your business' network. To access all of your services, all you need is a local network and an ADSL connection.

Naturally, it is immediately compatible with several Essilor tools, and there is no need for you to replace your equipment.

Designed around an open communication system, Essibox will gradually be interfaced with various management software and refraction or measurement equipment meeting international standards (please contact your local representative for details).

As your complete electronic assistant, Essibox blends in perfectly with the "architecture" of your organisation.

Flexible and upgradable, it enables you to gradually add new functions and services at your leisure and in line with your needs.

**ADAPTABILITY
+
FLEXIBILITY
+
UPGRADABILITY**



3. Essibox allows for the easy management of your jobs

Working in sync based around Essibox, your equipment will automatically transfer useful data (shape, prescription, height, PD, lens type, etc.), to the centring system and the edger, with no manual data entry required:

- The possibility to connect a centralised workshop serving several stores.
- The ability to directly prepare your choice of finishing on Essibox application: bevel position, groove width, chamfering....

- Visualisation of all of your jobs and their various states of progress.
- The possibility to organize or modify your priorities for the day.

Completely integrated within the shop's network, Essibox application is accessible from any computer workstation in your store.



CONSTANT VISIBILITY OF JOBS UNDERWAY


CONNECTION TO A CENTRALISED WORKSHOP SERVING SEVERAL STORES



4. With Essibox, create your own unique shapes and drilling patterns.

More than just a simple modification tool, Essibox also includes a shape and drilling creation module with which you can give free range to your imagination, bringing you access to new features reserved up until now for specialized companies.

Designed based on the underlying principles used for computer aided design software, this module stands out for its user friendliness and its new interface. The functions grouped together in the centre of the screen offer easy access.



SHOW OFF YOUR TALENTS TO BEST EFFECT!

PRODUCE CUSTOM SHAPES



5. Essibox manages your "shapes & drilling patterns" databases

Thanks to its unlimited storage capacity, Essibox enables you to create an archive of shapes and drilling models in complete freedom. The selection and identification of libraries are fast and simple thanks to:

- Multi-criteria searches (search by manufacturer, by collection or by model).
- Immediate visualisation of the shape and the related drilling pattern.

In addition to your own designs, you also have an integrated, updatable database for immediate use on Kappa c.t.D (which already includes 300 models from among the main brands).

If required, (and if the management software allows for this), you can also share your shapes to add to your customer files.



UNLIMITED
STORAGE
CAPACITY

INTEGRATED
DATABASE



6. Essibox optimises your machine usage

To keep your system running smoothly:

- Ongoing software updates for your digital system, for Essibox application and for the databases are available for download.
- Remote diagnosis service: warning messages available automatically or at your request for faster processing.
- Remote calibration (or guided calibration) for the equipment to optimise its performance.

- Custom service contracts: guaranteed intervention time, preventive maintenance visits, hotline, preferential terms for spare parts, etc...
- Availability of electronic documentation concerning your Essibox system and line.
- Telephone support for all additional explanations or specific help required.
- Secure access to the Professional websites.

NB: No information concerning your clients is transferred to Essilor.

Details of the custom service contracts are available from your local Essilor representative.

Coming soon...

MORE
NEW PARTERS

INNOVATING
SALES TOOLS