


OST 200 - OST 250



THE BENCHMARK FOR
HIGH COMFORT DAILY USE



OST 200 - OST 250

REFRACTION UNIT

Optometry





Standard chart projector holder





Optional chart projector holder



Optional drawer

	OST 200	OST 250		
ROTATING CHAIR			POWER SUPPLY	110/230V-50/60Hz
Motorised (elevation path 18cm)	●	●	DIMENSIONS AND WEIGHT	1.25 (L) x 2.1 (P) x 1.8 (H)m
Tiltable back 45°	●	●		240kg
Tiltable back 90°	-	●		
Lift-up elbow rests	●	●	INTEGRATED FUNCTIONS	Auxiliary electric socket Variable LED lighting Sensitive control panel Chart Projector plate (Dim : 21 x 30cm)
Footrests	●	●		
SLIDING TABLETOP			OPTIONS	Chart projector arm Storage drawer Side unit body covers: Red (RAL 3001), Blue (RAL 5005) or White (RAL 9002)
Dimensions: 42 x 82cm Max Weight: 50kg	2 instruments 90° Rotation	●		
Motorised height adjustment (elevation path 18cm)				
VERSIONS	Right or left handed version			

Products are sold separately. As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.



Instruments Department
64 bis, avenue Aubert
94306 Vincennes Cedex France
Tel: 33 (0)1 72 94 71 00
Fax: 33 (0)1 72 94 70 27
www.essilor-instruments.com

Your local distributor: