



OST 200

HIGH COMFORT REFRACTION UNIT
FOR DAILY USE



OST 200 WITH CONFORT PACK*

Tilting backrest 45°



Keypad




Chart projector holder with variable lighting



OST 200 CONFIGURATION:

Motorized chair (elevation path 18 cm)
Sliding table top 90°
Lift-up elbow rest
Footrest
Fixed chinrest
Right or left handed version
Dual voltage (110V / 220V) TBD

*CONFORT PACK includes:

- Tilting backrest 45°
- Keypad
- Chart projector holder with Variable lighting

OPTIONS:

Near vision lamp
LCD screen support
Chart projector holder
Storage drawer (60 x 350 x 70 mm)

Side unit body covers:

- Red (RAL 3001)
- Blue (RAL 5005)
- White (RAL 9016)

SPECIFICATIONS

Chair maximum load:

150 kg

Sliding tabletop:

Dimensions: 42 x 82 cm

Max. weight: 40 kg

Versions:

Right or left handed version

Power supply:




110/230 V- 50/60 Hz

Dimensions:

1.26 (L) x 2.14 (P) x 1.90 (H) m

Weight:

240 kg



ESSILOR INSTRUMENTS

81 boulevard Jean-Baptiste Oudry
94000 Créteil
France
Tél. 33(0)1 49 80 62 80
www.essilor-instruments.com

