

OST 250C

A COMPACT REFRACTION UNIT
TO FIT YOUR NEEDS




OST 250C A COMPACT REFRACTION UNIT TO FIT YOUR NEEDS

OST 250C CONFIGURATION

- Motorized chair (elevation path 18cm)
- Elevation table
- Sliding table top 90°
- Lift-up elbow rests
- Footrests
- Keypad
- Dual voltage (110V / 220V)
- Ambient light

The **OST 250C** is a compact refraction unit, the company's smallest yet. Its clever and harmonious design ensures patient comfort and an optimized eye exam, even in tight spaces.

- Conveniently sized at only 1 meter wide, **OST 250C** fits into small or narrow spaces.
- Its articulated manual/motorized arm adjusts easily.
- An adjustable, motorized seat accommodates all patients, whether children or adults.
- Available for left or right handed people to suit the needs of your particular space.



SPECIFICATIONS

Chair maximum load:
150kg

Sliding tabletop:
Dimensions: 42 x 82cm
Max. weight: 40kg

Versions:
Right or left handed version

Power supply:
100-110 and 200-230 V AC- 50/Hz


Dimensions & weight:
1.26 (L) x 2.14 (P) x 1.90 (H)m- 240kg

OPTION CHOSEN AT THE ORDER

Manual phoropter arm (OST 012)



Motorized phoropter arm



COLORS AVAILABLE FOR LATERAL COVERS

Red OST257  Blue OST255  White OST256 

CE Conformity marking

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.



ESSILOR INSTRUMENTS
81 boulevard Jean-Baptiste Oudry
94000 Créteil
France
Tel. : 33 (0)1 49 80 62 80
www.essilor-instruments.com

