



VISUAL COMFORT FOR THE PATIENT.
EASE AND FLEXIBILITY FOR THE ECP.



VISION-C 600

POLARIZED CHART SCREEN






24" POLARIZED LCD SCREEN

Multiple high resolution charts developed specifically for visual testing, meeting all optometry standards.

SLEEK LOOK & FEEL

Compliment your equipment with a modern and stylish design.

WIDE RANGE OF TESTS



VISUAL ACUITY:

Landolt C, Letters, Illiterate E, LogMar, Kids Optotypes

SPECIAL TESTS:

Binocular & fusion, Stereopsis, Strabismus, Aniseikonia, Red/Green, Phoria

OTHER:

Colour Vision, Random function, Contrast sensitivity, etc.

SPECIFICATIONS

Screen size:

24" LED

Luminosity:

Max 250 cd/m²

Resolution:

1920 x 1080

Projection distance:

3m-8m

Chart window size:

293.2 x 521.3

Charts:

Refer to Vision-R manual

Dimensions:

393 (L) x 612 (W) x 71 (D)

Weight:

7.9kg

Power supply:

12V - 5A (AV 100-120V/220-240V)

Power consumption:

Less than 60W

CE Conformity marking

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.

ERGONOMIC CONTROL

Fully integrated with the Essilor phoropter range for ease of use and automated refraction.




CUSTOMIZABLE

Quick and easy to customize the chart display.




MULTIPLE SET UP OPTIONS


Select your refraction distance between 3-8m with direct or mirror mode.



Screen Size



Linear Polarization



IR&IP Connectivity



28mm aesthetically



ESSILOR INTERNATIONAL

147 rue de Paris, 94220 Charenton-le-Pont - France
Tél. : +33 (0)1 49 80 62 80
www.essilor-instruments.com

